

New generation

# DrugRead®

## User-friendly interpretation of DrugWipe® tests



### Why use a Reader?

- Objective and documented DrugWipe® results
- Mobile and reliable drug screening
- Timer-controlled test run

### Advantages of DrugRead®

- Offers precise time keeping
- Shows a clear yes or no result
- Quality self test is integrated

### Site of operation

- Roadside control
- Checkpoints
- Police station

### Main Features:

Instrument:	Opto-electronical photometer for analysis and documentation
Housing:	Robust design and solid technology
Time control:	Handheld device with acoustical and optical alarm
User interface:	Easy to use touch screen with high brightness
Software / Language:	Customisable for each country
Battery performance:	approx. 300 measurements
Data storage:	64 MB for about 5.000 results
Interface:	Bluetooth interface for direct print out
Calibration:	Automatic unit self-test; additional calibration tests
Temperature range:	From +5 ° to +40 ° Celsius
Size:	21 cm x 10 cm x 12 cm
Weight:	~ 800 Gramm



we detect to protect